



SleepAngel distributor product catalogue
Consumer and hospitality range

SleepAngel – Where the Air is Pure

Sleep – rest and recovery – is crucial for our physical and mental performance and quality of life. At SleepAngel, we're committed to helping each individual get a comfortable and soothing sleep each and every night.

Our story started when our inventor David Woolfson's daughter got an infection from an unclean pillow. Seeing her suffer motivated him to research cross-infection. He realised that pillows were the most critical part because they are just little mattresses for your head where you are most vulnerable. After all, that's where you breathe, and it makes contact with eyes and ears, nose and mouth—all points of infection ingress.

Dr. Duncan Bain and David Woolfson invented PneumaPure Filter technology to stop infection transmission in healthcare setting. It is the world's first soft surface solution that combines absolute barrier protection (clinically proven to block pathogens*, allergens, microbes and mold) with user comfort.

SleepAngel bedding is scientifically proven, safe, sustainable, and supportive.

Product range includes pillows, mattresses, positioners, duvets, toppers and other soft surface products. SleepAngel for improved intimate bed environment hygiene at private homes and in hospitality industry.

Patent protected technology is independently tested and proven by Airmid Healthgroup.

*Tested pathogens include MRSA, C.Difficile, E.coli, Pseudomonas aeruginosa, Candida albicans, Aspergillus niger, Norovirus, Influenza Type A, Adenovirus 5, Coronavirus 229E (HCoV 229E)



PneumaPure Filter Technology Brief

Abstract

Design of bedding products (pillows, mattresses, positioners, duvets, top mattresses etc.) has not provided an adequate barrier to internal product contamination. Thread, needle perforations, air vents, zippers etc. allow ingress of air and liquid borne pathogens. Contaminated bedding acts both as a reservoir and potential vector for infection transmission and trigger for the sleeper.

Barrier bedding with PneumaPure™ filter technology

SleepAngel is the world's first bedding range to combine absolute barrier protection with airflow – for user comfort.

Combination of innovations – SafeWeld seams, PureShield fabric protection, mechanical PneumaPure filter – result in a barrier bedding product that blocks pathogens, bacteria, viruses, fungus and ingress of liquids, while remaining comfortable and supporting.



A white SleepAngel pillow is shown with three callout circles. One circle highlights the 'PneumaPure' filter technology on the top right. Another circle highlights the main body of the pillow. A third circle highlights a seam at the bottom right.

SleepAngel

Breathable mechanical filter that allows clean air flow into the pillow but blocks out pathogens and allergens.

Durable, soft touch and high performance textile that blocks out contaminants.*
Vapour permeable and waterproof.

Heat sealed seams provide complete barrier – no stitching holes.

* Tested pathogens include Coronavirus 229E (HcoV 229E), MRSA, E.Coli, Adenovirus 5, Pseudomonas aeruginosa, C.Difficile, Candida albicans, Aspergillus niger, Norovirus, Influenza Type A



Certified as Class I Medical devices in the EU
Registered Class I Medical Device with the FDA

SleepAngel®



Contents



6 HOME BEDDING RANGE

- 7 MICROFIBER PILLOWS
- 9 MEMORY FOAM PILLOWS
- 11 ORTHOPEDIC PILLOWS
- 13 TOPPERS
- 15 DUVETS
- 17 MATTRESSES

19 BABY & CHILD BEDDING RANGE

- 20 CHILDREN'S PILLOWS
- 22 CHILDREN'S DUVETS
- 24 CHILDREN'S TOPPERS
- 26 CHILDREN'S MATTRESSES

28 PERFORMANCE PILLOW RANGE

- 29 PERFORMANCE PILLOWS
- 32 MOVISTAR TEAM PERFORMANCE COLLECTION

34 TRAVEL PILLOW

- 35 TRAVEL PILLOW

37 HOSPITALITY BY SLEEPANGEL

HOME bedding range

Clean and comfortable recovery time



Microfiber pillows 1/2



Filling:	100% Polyester microfiber
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	1S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

SleepAngel®

Microfiber pillows 2/2



Size (cm)	Name	Filling weight (g)	Product code	Packaging per box*	MOQ (box)
50x60	Soft	500	3020109	5	1
50x60	Medium Firm	700	3020034	5	1
50x70	Soft	700	3020030	5	1
50x70	Medium Firm	900	3020014	5	1



*Carton dimensions 78x40x60 cm

Memory foam pillows 1/2



Filling:	100% Polyurethane foam 45 kg/m ³ (molded foam)
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	1S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Memory foam pillows 2/2



Direct foam size (cm) / Cover size (cm)	Name	Product code	Packaging per box*	MOQ (box)
42x56x11 / 50x60	Memory foam small	3020023	5	1
42x66x11 / 47x71	Memory foam medium	3020004	5	1

*Carton dimensions 78x40x60 cm

SleepAngel.

Orthopedic pillows 1/2



Filling:	100% Polyurethane foam 45 kg/m ³ (molded foam)
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	1S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Orthopedic pillows 2/2




Direct foam size (cm)	Name	Product code	Packaging per box*	MOQ (box)
50,5x33x6/10	Orthopedic Memory foam low	3020012	6	1
50,5x33x8/12	Orthopedic Memory foam medium	3020031	6	1
60,5x33x6/10	Orthopedic Memory foam low	3020027	5	1
60,5x33x8/12	Orthopedic Memory foam medium	3020029	5	1

*Carton dimensions 78x40x60 cm


SleepAngel®

Toppers 1/2



Filling:	100% Polyurethane foam, different foam types
Fabric:	Top Polyurethane coated nylon, 195 g/m ² , white; Bottom Antislip grey
Filter:	S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Toppers ^{2/2}



Size (cm)	Name	Product code	Packaging per pallet	MOQ (pallet)
80x200x5	Foam topper 80x200x5	3050100	5	1
90x200x5	Foam topper 90x200x5	3050101	5	1
160x200x5	Foam topper 160x200x5	3050104	5	1
80x200x5	Memory foam topper 80x200x5	3050110	5	1
90x200x5	Memory foam topper 90x200x5	3050111	5	1
160x200x5	Memory foam topper 160x200x5	3050114	5	1

SleepAngel®

Duvets 1/2



Filling:	100% polyester thermobonded wadding 200 g/m ²
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Duvets ^{2/2}



Size (cm)	Name	Product code	Packaging per box*	MOQ (box)
135x200	Duvet	4020001	6	1
150x200	Duvet	4020002	6	1
220x200	Duvet	4020025	4	1

*Carton dimensions 77x57x48 cm



Mattresses 1/2



Filling:	100% Polyurethane foam, different foam types
Fabric:	Top Polyurethane coated nylon, 195 g/m ² , white; Bottom Antislip grey
Filter:	S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.



Mattresses 2/2



Size (cm)	Name	Product code	Packaging per pallet	MOQ
140x200x25	Mattress	3080002	5	5 pcs
160x200x25	Mattress	3080004	5	5 pcs
200x200x25	Mattress	3080006	5	5 pcs

BABY & CHILD bedding range

Practical comfort for the smallest ones



SleepAngel.

Children's pillow 1/2



Filling:	100% Polyester microfiber
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	1S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Children's pillow 2/2



Size (cm)	Name	Filling weight (g)	Product code	Packaging per box*	MOQ (box)
40x60	3+ Year old children pillow	320	3020009	6	1

*Carton dimensions 78x40x60 cm

Children's duvet ^{1/2}



Filling:	100% polyester thermobonded wadding 200 g/m ²
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	1S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Children's duvet ^{2/2}



Size (cm)	Name	Product code	Packaging per box*	MOQ (box)
110x125	3+ Year old children duvet	4020024	10	1

*Carton dimensions 77x57x48 cm

Children's toppers ^{1/2}



Filling:	100% Polyurethane foam, different foam types
Fabric:	Polyurethane coated nylon, 195 g/m², white
Filter:	1S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Children's toppers ^{2/2}



Size (cm)	Name	Product code	Packaging per pallet	MOQ (pallet)
60x120x4	Nursery topper	3090021	16	1
70x140x4	Nursery topper	3090020	16	1

Children's mattresses 1/2



Filling:	100% Polyurethane foam, different foam types
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	1S filter 9,5x27 cm
Label:	Light grey print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Children's mattresses 2/2



Size (cm)	Name	Product code	Packaging per pallet	MOQ (pallet)
60x120x11	Nursery mattress	3090002	14	1
70x140x11	Nursery mattress	3090004	10	1

PERFORMANCE pillow range

Protecting your performance as you sleep



SleepAngel.

Performance pillows 1/2



Filling:	Microfiber, GelFlex Foam or Memory Foam
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	2 filters
Label:	Black print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Performance pillows ^{2/2}

Size (cm)	Name	Filling weight (g)	Product code	Packaging per box*	MOQ (box)
50x60	Microfiber Performance pillow 50x60 low	500	3020201	8	1
50x60	Microfiber Performance pillow 50x60 medium	700	3020202	8	1
50x70	Microfiber Performance pillow 50x70 low	700	3020203	8	1
50x70	Microfiber Performance pillow 50x70 medium	900	3020204	8	1
47x70x11	GelFlex Performance pillow 47x70 low		3020207	8	1
45x65x14	GelFlex Performance pillow 45x65 medium		3020208	8	1
47x71x11	Memory foam Performance pillow 47x71		3020209	6	1
50x60x11	Memory foam Performance pillow 50x60		3020210	8	1

*Carton dimensions 77x57x48 cm

SleepAngel is the official bedding partner for the Movistar Team.

The cooperation has inspired us to create Movistar Team custom design pillow collection and bring more importance to the role of rest and recovery in performance.

Custom design sets are available for minimum 100 set orders, please contact us for details and pricing.



Movistar Team Performance collection ^{1/2}



Filling:	Microfiber or GelFlex Foam
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	2 filters
Label:	Black print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

Movistar Team Performance collection ^{2/2}



Size (cm)	Name	Filling weight (g)	Product code	Packaging per box*	MOQ (box)
50x70	Microfiber Movistar Team 50x70 low	700	3020203 636 SET	8	1
50x70	Microfiber Movistar Team 50x70 medium	900	3020204 636 SET	8	1
47x70x11	GelFlex Movistart Team 47x70x11		3020207 636 SET	6	1

*Carton dimensions 77x57x48 cm

TRAVEL pillow

Clean place to rest your head wherever you are



SleepAngel.

Travel pillow 1/2



Filling:	GelFlex Foam
Fabric:	Polyurethane coated nylon, 195 g/m ² , white
Filter:	2 filters
Label:	Black print on white fabric
Cleaning:	Avoid abrasive scrubbing of the PneumaPure™ Filter. Do not dry clean. If soiled, wash under cold running water. For general hygiene, wipe down with a damp cloth when changing bed linen.
Safety Standards:	Cover flammability in accordance with BS 7175 source 5. UK filling qualities on request.

SleepAngel.

Travel pillow 2/2



Size (cm)	Name	Product code	Packaging per box*	MOQ (box)
28x47x11	Travel pillow	6020023	30	1

*Carton dimensions 77x57x48 cm

HOSPITALITY by SleepAngel

Safe, sustainable and supportive





Improved guest and staff safety with SleepAngel Barrier bedding – barrier for viruses, allergens, pathogens and liquids

Did you know that the pillow (and soft surfaces overall) is an often overlooked link in a chain of cross-infection

Traditional bedding requires add-on protecting covers. This increases your purchase & ongoing laundry costs


SleepAngel Barrier bedding for hospitality helps to address the Covid-19 challenges

SleepAngel PneumaPure™ technology, is globally patented (incl. EU, UK, North America, Japan, China, Russia, India, Middle East). Products are certified as Class I Medical Devices in the EU and registered Class I Medical Device with the FDA.

SleepAngel Hospitality Bedding range includes pillows, toppers, duvets and mattresses. Different fillings and sizes are available

	SleepAngel Bedding	Traditional Bedding
User comfort	Yes	Yes
Protected against Human Coronavirus	Yes	No
Asthma and Allergy friendly	Yes	N/A
Stays clean inside, barrier to Viruses, Allergens, Bacteria & Liquids	Yes	N/A
Easy to clean and maintain	Yes, just wipe clean	N/A
Periodic laundry wash	Not Required	Required
Warranty	2 Years	N/A


SleepAngel.



Homepage
www.SleepAngel.eu



VIDEO
SleepAngel - Sleep Where the Air is Pure
<https://youtu.be/malaNmVGSbw>



VIDEO
SleepAngel product technology
<https://youtu.be/GKyfgmQirrA>

SleepAngel.

Alejandro Valverde, Movistar Team:

"I am very pleased that Movistar Team partnered with SleepAngel to give us the best possible kit for recovery. SleepAngel Team specified a pillow for each of us that fits our personal comfort level. A comfortable and hypoallergenic pillow is a massive improvement for our rest and recovery time."

Stuart Hogg, Scottish rugby:

"I'm thoroughly enjoying having my SleepAngel pillow. It's become part of my travel pack and I take it everywhere with me. Great consistency in sleep, no matter where I am."

Juan Pablo Molinero, Movistar Team CMO:

"Sleep environment had been an outlier for some time with regard to controlling variables and protecting performance. We challenged SleepAngel to deliver on their promise of sleeping "Where the Air is Pure" and helping our athletes to recover better & reduce their risk of illness. We are delighted with the relationship and products so far and look forward to continuing our partnership. SleepAngel has become a must inside our protocols for racing."

